

VECTASHIELD® Hard♦Set™ Mounting Medium Compatibility

Based on in-house testing and data from published papers, Vectashield® Hard♦Set™ Mounting Medium has been used with the following fluorochromes.

Absence of a particular fluorochrome does not indicate incompatibility.

1,1'-dioctadecyl-3,3',3'-tetramethylindocarbocyanine	AMCA	hoechst dye
Alexa Fluor® 350	Cy™2	phycoerythrin
Alexa Fluor® 488	Cy™3	propidium iodide
Alexa Fluor® 546	Cy™5	rhodamine
Alexa Fluor® 555	DAPI	red fluorescent protein
Alexa Fluor® 568	DyLight® 488	SpectrumGreen™
Alexa Fluor® 594	DyLight® 549	SpectrumOrange™
Alexa Fluor® 633	DyLight® 594	SpectrumRed™
Alexa Fluor® 635	DyLight® 649	SYTOX® Green
Alexa Fluor® 647	fluorescein	tetramethylrhodamine
	green fluorescent protein	Texas Red®
		TO-PRO®-3

Trademarks are the properties of the respective companies.