

**Product Name:** Ponceau S Staining Solution

**Catalog No.:** 5225

**Batch No.:** 7

CAS Number: 6226-79-5

EC Number:

IUPAC Name: 3-Hydroxy-4-[2-[2-sulfo-4-[2-(4-sulfophenyl)diazenyl]phenyl]diazenyl]-2,7-naphthalenedisulfonic acid tetrasodium salt

**Description:**

Stain used to detect protein bands after Western blotting.


**Physical and Chemical Properties:**

Batch Molecular Formula: C<sub>22</sub>H<sub>12</sub>N<sub>4</sub>Na<sub>4</sub>O<sub>13</sub>S<sub>4</sub>

Batch Molecular Weight: 760.57

Physical Appearance: Red liquid

**Batch Molecular Structure:**



**Storage:** Store at RT

**Solubility & Usage Info:**

Soluble in 5% acetic acid/water (supplied as a 0.1% w/v solution)

**Stability and Solubility Advice:**

Some solutions can be difficult to obtain and can be encouraged by rapid stirring, sonication or gentle warming (in a 45-60°C water bath).

Information concerning product stability, particularly in solution, has rarely been reported and in most cases we can only offer a general guide. Our standard recommendations are:

**SOLIDS:** Provided storage is as stated on the product label and the vial is kept tightly sealed, the product can be stored for up to 6 months from date of receipt.

**SOLUTIONS:** We recommend that stock solutions, once prepared, are stored aliquoted in tightly sealed vials at -20°C or below and used within 1 month. Wherever possible solutions should be made up and used on the same day.

Caution - Not Fully Tested • Research Use Only • Not For Human or Veterinary Use

**bio-techne.com**

info@bio-techne.com

techsupport@bio-techne.com

**North America**

Tel: (800) 343 7475

**China**

info.cn@bio-techne.com

Tel: +86 (21) 52380373

**Europe Middle East Africa**

Tel: +44 (0)1235 529449

**Rest of World**

www.tocris.com/distributors

Tel: +1 612 379 2956